

Бутылки с узким и широким горлышком

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727) 34-54-704

Беларусь +(375) 257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: aqs@nt-rt.ru || сайт: <https://abdos.nt-rt.ru/>

Narrow Mouth Bottle, PP

- These bottles come with inner thread for superior sealing and better security.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.



Code	Capacity	Mouth Diameter	Height	No./case
PIII35	4 ml	12.2 mm	39.0 mm	72
PIII36	8 ml	17.5 mm	42.3 mm	72
PIII37	15 ml	17.5 mm	56.1 mm	72
PIII38	30 ml	17.5 mm	60.0 mm	72
PIII31	60 ml	17.5 mm	82.4 mm	72
PIII32	125 ml	21.0 mm	97.6 mm	72
PIII33	250 ml	21.0 mm	130.7 mm	72
PIII34	500 ml	24.6 mm	165.5 mm	48



Narrow Mouth Bottle, LDPE

- These bottles come with inner thread for superior sealing and better security.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.
- Suitable for storage of low-density liquids.



Code	Capacity	Mouth Diameter	Height	No./case
PIII45	4 ml	12.2 mm	39 mm	72
PIII46	8 ml	17.5 mm	42.3 mm	72
PIII47	15 ml	17.5 mm	56.1 mm	72
PIII48	30 ml	17.5 mm	60 mm	72
PIII41	60 ml	17.5 mm	82.4 mm	72
PIII42	125 ml	21 mm	97.6 mm	72
PIII43	250 ml	21 mm	130.7 mm	72
PIII44	500 ml	24.6 mm	165.5 mm	48



Narrow Mouth Bottle, HDPE

- These bottles come with inner thread for superior sealing and offers higher strength.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.
- Compatible with a wide range of laboratory reagents.



Code	Capacity	Mouth Diameter	Height	No./case
PIII55	4 ml	12.2 mm	39.0 mm	72
PIII56	8 ml	17.5 mm	42.3 mm	72
PIII57	15 ml	17.5 mm	56.1 mm	72
PIII58	30 ml	17.5 mm	60.0 mm	72
PIII51	60 ml	17.5 mm	82.4 mm	72
PIII52	125 ml	21.0 mm	97.6 mm	72
PIII53	250 ml	21.0 mm	130.7 mm	72
PIII54	500 ml	24.6 mm	165.5 mm	48



Amber Narrow Mouth Bottle, HDPE

- These bottles come with inner threads for superior sealing and offers higher strength.
- Opaque amber colour reduces UV light transmissions to protect light sensitive samples.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying guaranteed through the inclination of shoulders.
- Compatible with a wide range of laboratory reagents.



Code	Capacity	Mouth Diameter	Height	No./case
PIII55A	4 ml	12.2 mm	39.0 mm	72
PIII56A	8 ml	17.5 mm	42.3 mm	72
PIII57A	15 ml	17.5 mm	56.1 mm	72
PIII58A	30 ml	17.5 mm	60.0 mm	72
PIII51A	60 ml	17.5 mm	82.4 mm	72
PIII52A	125 ml	21.0 mm	97.6 mm	72
PIII53A	250 ml	21.0 mm	130.7 mm	72
PIII54A	500 ml	24.6 mm	165.5 mm	48



Narrow Mouth Bottle (Carboy Type), PP

- These bottles come with long inner threads for superior sealing and better security.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.



Code	Capacity	Mouth Diameter	Height	No./case
PIII27	2000 ml	24.9 mm	248.1 mm	6
PIII28	4000 ml	25.6 mm	333 mm	6
PIII29	8000 ml	38.3 mm	410 mm	6



Wide Mouth Bottle, PP

- These bottles come with inner threads for superior sealing and better security.
- Wide mouth of bottles allows for easy filling and emptying.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.



Code	Capacity	Mouth Diameter	Height	No./case
PIII60	30 ml	21.2 mm	58.9 mm	72
PIII61	60 ml	24.8 mm	81 mm	72
PIII62	125 ml	32 mm	94.5 mm	72
PIII63	250 ml	32.2 mm	126.7 mm	72
PIII64	500 ml	47.5 mm	163.2 mm	48



Wide Mouth Bottle, LDPE

- These bottles come with inner threads for superior sealing and better security.
- Wide mouth of bottles allows for easy filling and emptying.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.
- Suitable for storage of low-density liquids.



Code	Capacity	Mouth Diameter	Height	No./case
PII170	30 ml	21.2 mm	58.9 mm	72
PII171	60 ml	24.8 mm	81 mm	72
PII172	125 ml	32 mm	94.5 mm	72
PII173	250 ml	32.2 mm	126.7 mm	72
PII174	500 ml	47.5 mm	163.2 mm	48



Wide Mouth Bottle, HDPE

- These bottles come with inner threads for superior sealing and better security.
- Wide mouth of bottles allows for easy filling and emptying.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.
- Compatible with a wide range of laboratory reagents.



Code	Capacity	Mouth Diameter	Height	No./case
PII180	30 ml	21.2 mm	58.9 mm	72
PII181	60 ml	24.8 mm	81 mm	72
PII182	125 ml	32 mm	94.5 mm	72
PII183	250 ml	32.2 mm	126.7 mm	72
PII184	500 ml	47.5 mm	163.2 mm	48



Amber Wide Mouth Bottle, HDPE

- These bottles come with inner thread for superior sealing and better security.
- Opaque amber colour reduces UV light transmissions to protect light sensitive samples.
- Wide mouth of bottles allows for easy filling and emptying.
- These bottles are used for packing, shipping and storage purposes.
- Superior resistance to chemicals.
- Complete emptying is guaranteed through the inclination of shoulders.
- Compatible with a wide range of laboratory reagents.



Code	Capacity	Mouth Diameter	Height	No./case
PII180A	30 ml	21.2 mm	58.9 mm	72
PII181A	60 ml	24.8 mm	81 mm	72
PII182A	125 ml	32 mm	94.5 mm	72
PII183A	250 ml	32.2 mm	126.7 mm	72
PII184A	500 ml	47.5 mm	163.2 mm	48



Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727) 34-54-704

Беларусь +(375) 257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

По вопросам продаж и поддержки обращайтесь:

эл.почта: aqs@nt-rt.ru || сайт: <https://abdos.nt-rt.ru/>